


MERYart is a line of modern premium bed head units where medical functionality and high design esthetics come together. Available with Corian® or tempered glass fronts, MERYart accommodates medical gases, electrical and IT systems. Customisable graphics are available.



TECHNICAL SPECIFICATION

- Body made of anodized or powder coated aluminium profile mounted directly on the wall
- Decorative front panel made of Corian®

Size:

- Length (standard): single-unit: 1500 mm, double unit: 3200 mm, triple unit: 4800 mm
- Depth: 77 mm
- Height: 370 mm

MEDICAL GAS OUTLETS AND VACUUM

- Available standards: SS / DIN / BS / AFNOR / others
- Medical gases location on the BHU: at the door side

LIGHTING

- Lighting types: LED / T5 / 2G11
- Lighting configuration: general/room, reading, night
- Lighting is switched on/off by switch/button on the BHU and/or the wall and/or patient's handset

CHARACTERISTICS

- Medical gas outlets, electrical sockets and lighting switches are located on the front side of the unit, similar to MERYlight type
- Graphics selection customized depending on customers demand

STANDARD CONFIGURATION IN HOSPITAL WARDS



- Medical gas outlets: 1× O₂, 1× VAC
- 230 V electrical sockets: 3 pcs (2 circuits)
- Equipotential pins: 1 pc
- IT socket RJ-45, Cat. 6: 1 pc
- Reading lighting 1× 12 W LED, 2200 lm; switch on the unit
- Night lighting 1× 4 W LED, 300 lm; switch on the unit
- NCS cut-out

ELECTRICAL AND IT SOCKETS

- Electrical sockets: Type E (FR) / Type F (Schuko) / Type G (BS) / other
- IT sockets: RJ-11 / RJ-12 / RJ-45, 6a cat. or higher
- Multimedia sockets: USB, HDMI, other
- NCS depending on customers' demand

VERSION A

E – Electrical system inlet, 180 x 50 mm
G – Gas system inlet, 200 x 50 mm



VERSION B

(mirror reflection)

