

MERYfun is a bed head unit designed for children's wards. It combines medical functionality with a friendly form. Thanks to decorative Corian® fronts with colourful graphics, it helps to create a cheerful and welcoming space for the youngest patients. Integration of lighting, all media and a nurse call system ensure that this medical supply unit is fully functional.



TECHNICAL SPECIFICATION

- Body made of anodized aluminium profile mounted directly on the wall
- The decorative front panel made of Corian®.

Size:

- Length (standard): single-unit: 1500 mm
- Depth: up to 100 mm, depending on the decorative panel design
- Height: 370 mm

MEDICAL GAS OUTLETS AND VACUUM

- Available standards: SS / DIN / BS / AFNOR / others
- Medical gases location on the BHU: at the door side

LIGHTING

- Lighting types: LED / 2G11
- Lighting configuration: general/room, reading, night
- Lighting is switched on/off by switch/button on the BHU and/or the wall and/or patient's handset

CHARACTERISTICS

- Medical gas outlets, electrical sockets and lighting switches are located on the front side of the unit, similar to MERYlight type
- Decorative front panel design selected from the available design library

STANDARD CONFIGURATION IN HOSPITAL WARDS

- Medical gas outlets: 1× O₂, 1× VAC
- 230 V power sockets: 3 pcs (2 circuits)
- Equipotential pins: 1 pc
- IT socket RJ-45, cat. 6: 1 pc
- Reading lighting 1× 12 W LED, 2200 lm; switch on the unit
- Night lighting 1× 4 W LED, 300 lm; switch on the unit
- NCS cut-out

ELECTRICAL AND IT SOCKETS

- Electrical sockets: Type E (FR) / Type F (Schuko) / Type G (BS) / other
- IT sockets: RJ-11 / RJ-12 / RJ-45, 6a cat. or higher
- Multimedia sockets: USB, HDMI, other
- NCS depending on customers' demand

VERSION A

E – Electrical system inlet, 180 x 50 mm
G – Gas system inlet, 200 x 50 mm

