

**LISSAclassic** ceiling pendant is designed to distribute gas and electrical media and to provide ergonomic access to equipment necessary for the work of medical personnel.

This pendant is used in:

- Operating rooms
- Intensive Care units, Critical Care units
- Emergency units
- Treatment rooms

Depending on configuration LISSAclassic can function as:

- **anaesthesia pendant**
- **surgical pendant** (for various specialisations)
- **ICU pendant**

### ARMS TECHNICAL SPECIFICATION

- Single/double swivel arm
- Single arm lengths: 600, 800, 1000, 1200 mm
- Total arms range: 1200, 1400, 1600, 1800, 2000, 2200, 2400 mm
- Possibility of using non-standard arms ranges
- Rotation – maximum rotation range: 350°, can be reduced by adjusting rotation limiters
- Brake: pneumatic, electromagnetic or friction
- Arms, ceiling covers, service head covers and accessories powder coated in RAL 9002
- Arms end-caps, drop tubes, bottom covers and blind covers powder coated in RAL 9022

### SERVICE HEAD TECHNICAL SPECIFICATION

- Horizontal service heads' width: 600, 700, 800, 900, 1000 mm
- Vertical 3, 4 or 6-plane service heads length: 600, 800, 1000, 1200, 1500 mm
- Non-standard service head size possible
- Service head's body made of aluminium profiles
- Non-standard RAL for service head covers possible

### MEDICAL GAS AND VACUUM OUTLETS

- Standard: SS / DIN / BS / AFNOR / other
  - Type: O<sub>2</sub>, AIR, N<sub>2</sub>O, CO<sub>2</sub>, VAC
- Gas outlets quantities upon to individual configuration
- AGSS outlet

### ELECTRICAL AND IT SOCKETS


- Type: Type E (FR) / Type F (Schuko) / Type G (BS) / other
- IT sockets: RJ-11 / RJ-12 / RJ-45, 6a cat. or higher
- Multimedia sockets: USB, HDMI and other
- Nurse call system upon to user's needs

### OPTIONAL LIGHTING

- in upper part of ceiling arm
- in side cover of service head
- in the bottom of service head
- mobile (individual lamp), fixed on the rail

### OPTIONAL EQUIPMENT

- Vertical poles for accessories
- Medical rails made of acid-resistant steel
- Additional accessories:
  - different size and load capacity shelves
  - drawers
  - IV poles
  - monitor hangers / arms
  - catheters and medical disposables baskets
  - cables and hoses organizers



#### Vertical head sections

