

**LUNAcompact** ceiling pendant is intended to gas and electricity distribution and provides ergonomic access to equipment necessary for the work of medical personnel. This minimalistic solution has been designed to provide access to media within a strictly defined space.

This pendant is used in:

- **Intensive Care Units**
- **Treatment and dressing rooms**
- **Hospital rooms with lighter equipment requirements**
- **Auxiliary and preparation rooms**

Depending on configuration **LUNAcompact** can function as:

- **anaesthesia pendant**
- **treatment pendant**
- **ICU pendant**

Due to its construction, LUNAcompact is recommended for use as **ICU pendant** in less demanding hospital wards where pendants mobility is not required.

## SERVICE HEAD TECHNICAL SPECIFICATION

- Horizontal service heads width: 600, 700, 800, 900, 1000 mm
- Vertical 3, 4 or 6-plane service heads length: 600, 800, 1000, 1200, 1500 mm
- Non-standard size service heads possible
- Service head's body made of aluminium profiles
- Ceiling covers, service head covers and accessories powder coated in RAL 9002
- Bearings end-caps, drop tubes, bottom covers and blind covers powder coated in RAL 9022
- Non-standard RAL for service head covers possible

## MEDICAL GAS AND VACUUM OUTLETS

- Standards: SS / DIN / BS / AFNOR / other
  - Type: O<sub>2</sub>, AIR, N<sub>2</sub>O, CO<sub>2</sub>, VAC
- Gas outlets quantities upon to individual configuration
- AGSS outlet

## ELECTRICAL AND IT SOCKETS

- Electrical sockets: Type E (FR) / Type F (Schuko) / Type G (BS) / other
- IT sockets: RJ-11 / RJ-12 / RJ-45, 6a cat. or higher
- Multimedia sockets: USB, HDMI and other
- Nurse call system upon to user's needs

## OPTIONAL LIGHTING

- in side cover of service head
- in the bottom of service head
- mobile (individual lamp), fixed on the rail

## OPTIONAL EQUIPMENT

- Vertical poles for accessories
- Medical rails made of acid-resistant steel
- Different size and load capacity shelves
- Drawers
- IV poles
- Monitor hangers / arms
- Catheters and medical disposables baskets
- Cables and hoses organizers



### Vertical head sections

