

MERYcare Classic is an advanced unit coming in horizontal and vertical versions and designed to be used in Intensive Care Units (ICU/NICU). It combines high quality, functionality and ergonomics to support medical staff in work requiring fast action, precision and reliability.

TECHNICAL SPECIFICATION

- Body made of anodized or powder-coated aluminium profile mounted directly on the wall
- Front cover powder-coated in any RAL color
- Number of inlets (gas / electrical systems) depending on configuration and fittings

Size:

- Length (standard): Single unit – 2000 mm,
- Non-standard lengths available
- Depth: 110 mm
- Height: 240 / 270 / 370 / 400 mm – depending on configuration

MEDICAL GAS OUTLETS AND VACUUM

- Available standards: SS / DIN / BS / AFNOR / other
- Medical gases location on the BHU: at the door side

LIGHTING

- Lighting types: LED / 2G11
- Lighting configuration: general/room, reading, night
- Lighting is switched on/off by switch/button on the BHU and/or the wall and/or patient's handset

CHARACTERISTICS

- Electrical and IT sockets fixed at an angle of 25° to the floor

STANDARD CONFIGURATION IN HOSPITAL WARDS


- Single unit 2000 mm: recommended
- Medical gas outlets: 2x O₂, 2x AIR, 2x VAC
- 230 V electrical socket: 6 pcs. (2 circuits)
- Equipotential pins: 2 pcs
- IT sockets: RJ-45, 6 cat., 2 pcs
- General lighting: 1x 24W LED, 4400 lm, Ra>90: switch on the unit
- Reading lighting 1x 12W LED, 2200 lm, Ra>90: switch on the unit
- Night lighting 1x 4W LED, 300 lm: switch on the unit
- NCS cit-out
- Medical rail fixed on the unit: 2 pcs

ELECTRICAL AND IT SOCKETS

- Electrical sockets: Type E (FR) / Type F (Schuko) / Type G (BS) / other
- IT sockets: RJ-11 / RJ-12 / RJ-45, 6a cat. or higher
- Multimedia sockets: USB, HDMI, other
- NCS depending on customers' demand


VERSION A

E – Electrical system inlet 180 x 50 mm
G – Gas system inlet, 200 x 50 mm



VERSION B

(mirror reflection)



OPTIONAL EQUIPMENT

- Medical rails made of acid-resistant steel, fixed on the unit and/or wall
- Additional accessories: shelves (including shelves with drawers) of various sizes and load capacities, examinations lamps, amonitors holders, rods for infusion pumps, baskets for disposables

Medical device, Class IIb in accordance with Regulation (EU) 2017/745.