

**in**med

# LISSAclassic



SCAN TO SEE LISSAclassic IN ACTION



# LISSAclassic



## LISSA Ceiling Pendants – The Future Begins Today

- Quality • Reliability • Strong market position
- New generation = a new era of ceiling pendants
- Engineering excellence inspired by the future



## NEW Premium LISSAclassic Ceiling Pendant

- Rail + ergonomic handles: fast accessories mounting & smooth movement
- Smoothly integrated display: time + BMS alarms + FM + Bluetooth audio
- Environmental monitoring: noise level + TVOC air quality
- Multi-axis swivel arm + rotating shelf
- Retractable section: medical gases + power + IT + cable storage
- Dedicated configurations for OR / ICU / NICU
- Smooth movement even when fully loaded

Medical device, Class IIb in accordance with Regulation (EU) 2017/745.



INMED S.A. Kwiatowa 32A 55-330 Krępiec, POLAND  
+48 71 398 02 28; [www.inmed.pl/en](http://www.inmed.pl/en)

This document does not constitute a commercial offer.

Leaflet version: EN/2026/LISSAclassic/1