

Bed Head Unit



SCAN TO SEE MERYfun MBAND VIDEO

MERYfun MBAND BHU – PIANO

Designed for the youngest patients, it combines medical functionality with play and therapy. It ensures full safety of use – electrical sockets are hidden when not in use, and access to the unit's functions is protected by electronic keys. The solid Corian® housing by DuPont™ guarantees durability and hygiene. **MERYfun PIANO** is more than medical supply unit – it is an interactive device that plays music!

Key features

- Modern design
- Coherent esthetics that goes beyond functionality of this type of units
- Flush-mounted electrical socket module with the option to hide sockets
- Interactive piano keyboard allowing play like a real instrument
- Resonance system using the front housing as a speaker
- Color samples according to Corian color range

Technical specification

- Dimensions: 918 x 1008 x 77 mm
- Material: Corian Deep Nocture / Glacier White, aluminium mounting frame
- Hidden electrical sockets: 45x45 system
- IT sockets: RJ-11 / RJ-12 / RJ-45 (Cat. 6a+)
- Multimedia: USB

Optional accessories: Wall-mounted medical rails made of acid-resistant steel, shelves of various sizes & drawers fixed to medical rails, examination lamps, monitor articulated arms, IV poles, baskets for medical disposables

Optional designs: drums, guitar.



Example configuration:

- 1x O2, 1x VAC, 3x 230 V, USB A+C, 1 module for user's own installation
- Medical gas outlet systems: SS / DIN / BS / AFNOR
- Optional lighting: LED Ambient Light 500 lm



Medical device, Class IIb in accordance with Regulation (EU) 2017/745.



INMED S.A. Kwiatowa 32A 55-330 Krępiec, POLAND
+48 71 398 02 28; www.inmed.pl/en

This document does not constitute a commercial offer.
Leaflet version: EN/2026/MERYfun_MBAND/1